

BeamOn *2/3" CCD Beam Profiler*



Main Features

- Real time beam size and Gaussian fit
- 2D/3D plots of beam in real time
- Software controlled electronic shutter & gain
- Video with playback, snapshot files
- Data exporting to another computer via RS232/TCP-IP
- Data logging with detailed statistics
- Automatic Pass/Fail analysis report
- ActiveX package to control software from your application

Specifications

- Spectral range: 350-1100nm
- Camera type: Monochrome interline transfer CCD 2/3" format
- Pixel size: 11.6 μm X 11.2 μm
- Measuring area: 8.5mm X 6.35mm
- Frame rate: 25 Hz
- External trigger
- Compact in size: 6mm (L) X 4.5mm (W) X 4mm (H)
- C-Mount aperture
- Operating systems: Windows XP/2000/7 32 & 64 bit



DUMA OPTRONICS LTD.

1st Hazait St., P.O.Box 3370 Neshar 20306, Israel Tel:972-4-8200577 Fax: 972-4-8204190
 Website:<http://www.duma.co.il> E-mail:sales@duma.co.il