





MXC connector for CCD camera


Green led noting system's readiness



Optional power connector 5V DC

1.8 m USB cable



				DIMENSIONS IN MILLIMETERS (mm)		scale 1:1	
				 <b>DUMA OPTRONICS LTD.</b>			
				TITLE: <b>USB attachment for SPOTON CCD</b>			
				Do not scale drawing			
				DWG NO. <b>USB_A_01</b>			
				REVISION:			
	NAME	SIGNATURE	DATE				
DRAWN	Yair		28.2.2012				
CHKD	Oren		28.2.2012				
APPVD	Judith		28.2.2012				