

# Adjusting X-Y Stage for EAC Line



- Optional Adjustment Stage, driven by two micrometers, to be mounted on the Autocollimator's basis.
- Compatible with the following models:
  - EAC-1012-19
  - EAC-1012-L19
  - EAC-1012-19-FO
  - EAC-1012-19-XR

**Dimensions:** 156 x 110 x 47 [mm]

Out-of-plane tilting range:  $\pm 2^\circ$

In-plane rotation range:  $\pm 2.5^\circ$



## Ordering Information

**EAC-1012-19-STG:** Optional Accessory for XY Table adjustment

**DUMA OPTRONICS LTD.**

Website: <http://www.dumaoptronics.com>

E-mail: [sales@duma.co.il](mailto:sales@duma.co.il)

July 2022

