

Pan/Tilt Adjusting Stage for EAC Line



- Optional Adjustment Stage, driven by two micrometers, to be mounted on the Autocollimator's basis.
- Compatible with the following models:
 - EAC-1012-19
 - EAC-1012-L19
 - EAC-1012-19-FO
 - EAC-1012-19-XR

Dimensions: 156 x 110 x 47 [mm]

Out-of-plane tilting range: $\pm 2^\circ$

In-plane rotation range: $\pm 2.5^\circ$



Ordering Information

EAC-1012-19-STG: Optional Accessory - Pan/Tilt Stage.

DUMA OPTRONICS LTD.

