


Flexible Filter Wheel



- Optional Flexible Filter Wheel, matching the Autocollimator's output aperture with dedicated Neutral Density Filters.
- Optional polarization filters are available.
- Has 5 mountable positions.
- Directly powered & controlled via USB through a dedicated GUI.
- Compatible with the following models:
 - EAC-1012-19
 - EAC-1012-L19
 - EAC-1012-19-FO
 - EAC-1012-19-XR
 - Beam Profiling Line-Up

Ordering Information

EAC-1012-19-FW5: Optional Accessory providing a 36 [mm] clear aperture filter wheel for ND Filters, Polarization Filters and others.

DUMA OPTRONICS LTD.

Website: <http://www.dumaoptronics.com>

E-mail: sales@duma.co.il

July 2023

