



# BeamAnalyzer Rotation Mount



The rotation mount is designed to attach to the M4 holes at the BeamAnalyzer head and provide precise rotation of the head about the optical axis.

It is sometimes necessary to rotate the BeamAnalyzer head about the optical axis. Each head is provided with three M4 threaded mounting holes on its rear plate for this purpose.

Dimensions (mm):        80 x 62 x 36

## Ordering Information

**BA-Mount:** A mount enabling head rotation about the optical axis.

**DUMA OPTRONICS LTD.**

Website: <http://www.dumaoptronics.com>

E-mail: [sales@duma.co.il](mailto:sales@duma.co.il)

October 2023

