


Flexible Filter Wheel



- Optional Flexible Filter Wheel, matching the Autocollimator's output aperture with dedicated Neutral Density Filters.
- Optional polarization filters are available.
- Has 5 mountable positions.
- Directly powered & controlled via USB through a dedicated GUI.
- Compatible with the following models:



- EAC-1012-19
- EAC-1012-L19
- EAC-1012-19-FO
- EAC-1012-19-FO-E/T
- EAC-1012-19-XR
- Beam Profiling Line-Up

Arrangement example of a filter wheel along with the autofocusing electronic autocollimator

Ordering Information

EAC-1012-19-FW5: Optional Accessory providing a 36 [mm] clear aperture filter wheel for ND Filters, Polarization Filters and others.

Filters included: ND2, ND4, Polarizing Filter (0°), Polarizing Filter (90°), one unpopulated.

DUMA OPTRONICS LTD.

Website: <http://www.dumaoptronics.com>

E-mail: sales@duma.co.il

January 2024

